

SALVO+GEL

INSTANT HAND SANITIZER



Germs remain active on hard surfaces for hours, and sometimes even days. These germs are passed on from one person to another when we handle everyday things like telephones, computer keyboards and mice, door handles, hand rails, money, etc. Experts acknowledge that approximately 80% of all common infections are spread by the hands this way.

SALVO-GEL is an important weapon in the ever more serious battle against pathogenic bacteria.

SALVO-GEL has been developed to meet the requirements of the food industry for a completely effective, yet quick-acting and nontainting, food-safe hand steriliser. Combining a high level of alcohol gives proven bactericidal effectiveness.

Added glycerine leaves hands soft, supple, non-sticky and moisturised - does not dry out the skin and leaves no after-odour.

These advanced properties make SALVO-GEL equally suitable for many medical and care service applications - in fact, any application where there is a risk of cross-infection from hand contact.

Independent laboratory tests prove 99.9999% kill of gram positive and gram negative bacteria.

Complies with EU Approval EN1276 for food use and is effective against *Pseudomonas Aeruginosa*, *Escherichia Coli* ATCC 10536, *Staphylococcus Aureus* ATCC 6538, and *Enterococcus hirae* after a contact time of 5 minutes.

Complies with EU Approval EN13727 for effective use against MRSA, *Salmonella Typhimurium*, *Listeria Monocytogenes* in medical areas.

Apply a small quantity to clean, dry hands and rub well in, paying particular attention to fingernails and creases in the skin. Within 15 to 20 seconds, SALVO-GEL will have evaporated, leaving the skin soft and dry. It may also be used in the same way to sterilise surgical gloves.

COSHH Data: An EC safety data sheet (MSDS) providing additional information is available on request.

Packaging Sizes:

- 24 x 80ml pump or flip-top
- 12 x 250ml pump or flip-top
- 10 x 500ml pump bottles
- 6 x 1 litre pump bottles
- 6 x 800ml pouches for use in 800ml pouch dispenser (illustrated)
- 6 x 2000ml pouches for use in 2 litre pouch dispenser
- 4 x 5 litre bulk-fill jerrycans

SALVO-GEL IS USED IN:

- Clinics
- Dental surgeries
- Doctors' surgeries
- Food preparation areas
- Hospitals
- Kindergartens
- Nursing homes
- Schools
- Colleges
- Butchers' shops
- Veterinary surgeries
- Local Councils
- Supermarkets
- Bakeries



Excalibur Hygiene Limited

Units 4/5 Dadsford Bridge Industrial Estate, Plant Street, Wordsley, Nr. Stourbridge, West Midlands DY8 5SY

Tel: 01384 671505 | Fax: 01384 671256

E-mail: sales@excalibur-hygiene.co.uk